

# HORIZONS INTERNATIONAL BIOPSY FORCEPS



P.O Box 213-3006 Zona Franca Metro, Barreal de Heredia, Costa Rica  
Phone (506) 2261-9114 • Fax (506) 2261-9115  
[www.horizonscorp.com](http://www.horizonscorp.com) • [info@horizonscorp.com](mailto:info@horizonscorp.com)



## BIOPSY FORCEPS

Horizon International biopsy forceps sever the toughest tissue and obtain the largest possible specimens with the least amount of effort. All moving parts are precision machined for ease of operation. Even the unique Teflon® inner tube is engineered for smooth, effortless action no matter how many twists and turns required to extract biopsy sample.

Choose from a complete range of forceps to fit all your gastroenterology and pulmonary needs, with features that include:

- Surgical steel jaws maintain a good cutting edge procedure after procedure.
- Forceps can be sterilized using ETO, autoclaves or soak method.
- Sizes to fit all your biopsy needs.
- Handle designed for tactile feel, comfort and ease of use.
- Clean, quick, economical and efficient.
- Manufactured and assembled in the U.S.A.
- Available with Teflon coated sheath.

P.O Box 213-3006 Zona Franca Metro,  
Barreal de Heredia, Costa Rica  
Phone (506) 2261-9114 • Fax (506) 2261-  
9115  
[www.horizonscorp.com](http://www.horizonscorp.com) • [info@horizonscorp.com](mailto:info@horizonscorp.com)



## BIOPSY FORCEPS

### Ordering Information

Product Application	Product Number	Jaw Description	Recommended Minimum Channel Dia.	Working Length
Pediatric Bronchoscope		Oval	1.2 mm *	115cm
		Alligator	1.2 mm *	115cm
Bronchoscope		Oval	2.0 mm **	100cm
		Oval-Spike	2.0 mm **	100cm
		Alligator	2.0 mm **	100cm
		Alligator-Spike	2.0 mm **	100cm
Therapeutic Bronchoscope Sigmoidoscope		Oval	2.5 mm	100cm
		Oval-Spike	2.5 mm	100cm
		Alligator	2.5 mm	100cm
Pediatric Gastroscope		Alligator-Spike	2.5 mm	100cm
		Oval	2.0 mm **	160cm
		Oval-Spike	2.0 mm **	160cm
		Alligator	2.0 mm **	160cm
Gastroscope		Alligator-Spike	2.0 mm **	160cm
		Oval	2.5 mm	240cm
		Oval-Spike	2.5 mm	240cm
		Alligator	2.5 mm	240cm
Colonoscope		Alligator-Spike	2.5 mm	240cm
		Oval	2.5 mm	240cm
		Oval-Spike	2.5 mm	240cm
		Alligator	2.5 mm	240cm
		Alligator-Spike	2.5 mm	240cm
		Oval	3.4 mm	240cm
		Oval-Spike	3.4 mm	240cm
	Alligator	3.4 mm	240cm	
	Alligator-Spike	3.4 mm	240cm	

**Optional Teflon® Exterior Coating.** Optional exterior Teflon coating offers the ultimate in endoscope channel protection, which greatly reduces downtime and repair cost. Passes easily through the instrument channel. To order Teflon coating, follow the product number with a “T” (e.g., 18100-1T). Call for more information.

\*Not available with Teflon coating.

\*\*Minimum channel diameter of 2.2mm is required for Teflon-coated forceps

P.O Box 213-3006 Zona Franca Metro,  
Barreal de Heredia, Costa Rica  
Phone (506) 2261-9114  
Fax (506) 2261-9115  
www.horizoncorp.com  
info@horizoncorp.com