

# COBRA™

The Power of Two in One  
Rotatable Snare / Retrieval Device



360° Rotation


Horizons International Corp.

*We Listen, Innovate, and Create to save lives.*

# COBRA™

## The Power of Two in One

### Rotatable Snare / Retrieval Device



FEATURES	BENEFITS
The Power of Two in One	Unique Patented design, half snare and half Retrieval Net, to cut and collect polyps in one single pass.
360° Rotation (1:1 ratio in all directions)	Designed to facilitate cutting, capturing and removing of polyps, in one single pass.
Deep, Three - Dimensional Net	Engineered to capture and remove multiple polyps, in a single pass.
Durability: Interwoven Net Design	Net is attached to the outer wire using an interwoven design to help maintain the integrity of the net.
Short throw handle	Ergonomically designed for user comfort.

#### PRODUCT SPECIFICATIONS

### COBRA™ Rotatable Snare / Retrieval Device

PRODUCT CODE	CATHETER DIAMETER MM	SIZE MM	LENGTH CM	TYPE	PACKAGING
HIC-SODT-13-5	2.5	13	230	OVAL	BOX 5
HIC-SODT-13-10	2.5	13	230	OVAL	BOX 10
HIC-SODT-24-5	2.5	25	230	OVAL	BOX 5
HIC-SODT-24-10	2.5	25	230	OVAL	BOX 10
HIC-SODT-25-5	2.5	32	230	OVAL	BOX 5
HIC-SODT-25-10	2.5	32	230	OVAL	BOX 10
HIC-SHDT-13-5	2.5	13	230	HEXAGONAL	BOX 5
HIC-SHDT-13-10	2.5	13	230	HEXAGONAL	BOX 10
HIC-SHDT-24-5	2.5	25	230	HEXAGONAL	BOX 5
HIC-SHDT-24-10	2.5	25	230	HEXAGONAL	BOX 10
HIC-SHDT-25-5	2.5	32	230	HEXAGONAL	BOX 5
HIC-SHDT-25-10	2.5	32	230	HEXAGONAL	BOX 10

COBRA™ SNARE is a Trademark of Horizons International Corp

Protected by one or more of the following US Patents: 9,101,342; 8,858,567; 10,966,748 and Pending Patents.

Caution: Federal (USA) law restricts this device to be sold by, or on the order of a physician.

Distributed by:



*We Listen, Innovate, and Create to save lives.*

**Horizons International Corp.**

P.O. Box 213-3006 Zona Franca Metro,  
Barreal de Heredia, Costa Rica

Phone: (506) 2261-9114/ Fax: (506) 2261-9115

**FOR ORDERS IN THE U.S. PLEASE CALL TOLL FREE 866-245-7697**